

VNMR 6.1C 1D NMR Quick Guide (Nightmare/NMR500)

1. To Begin:

- Logon to computer with your username and password
- VNMR Software should launch automatically, if not, Right Click (RC) on desktop, if you don't see xterm, RC on Program, then Left Click (LC) on shell tool. Then type *vnmr* on the UNIX prompt.
- Type *e* <rtn>. Place sample in the spinner. Gauge properly. Place on top of magnet. Type *i* <rtn>.
- Type *reshim* <rtn>. Enter solvent name, for e.g. *CDCl3* <rtn>.
- Click *Acqi* button. A new window will pop up. If you don't see the *Acqi* button, you can type *acqi* <rtn>.

2. Lock and Shimming:

- Check spinning, it should be off (set to 0).
- Click on **LOCK** panel, Click **LOCK off**.
- Increase lockpower and gain.
- **Lockphase** Nightmare=144; NMR500=310
- Move Z0 slider slowly until one 'beat' is visible. Reduce the lockpower and gain if lock level went over 100. Keep adjusting Z0, lockpower, and gain until the lock level maximized at around 70.
- Click **LOCK on**. Then click **SHIM**.
- Adjust Z1 with **± 16 button** to maximize lock level. Repeat on Z2 (Z2C for NMR500 and using **±4 button**). Adjust X1 follow by adjustment on Y1. Continue to alternate Z1, Z2, X1, and Y1 until no improvement.

3. Shimming:

- Click **CLOSE** to exit the lock/shimming window.
- Type *proton* <rtn> to load a proton experiment. Also enter the appropriate solvent name: *CDCl3* <rtn>.
- Type *nt=1 ss=0 ga* <rtn>
- When complete, type *f full aph* <rtn>, expand around solvent peak or suitable well-resolved singlet. Type *vsadj* <rtn>.
- Is this peak well shimmed (is it narrow or symmetric)? If yes, proceed to acquisition. If not, click *Acqi* then *SHIM*.
- Adjust appropriate shim (e.g. Z1 for symmetric broadening or Z2 for asymmetric peak shape).

4. Acquiring a Spectrum:

- Type *proton* for 1D 1H experiment and type *carbon* for 1D carbon. Also enter the appropriate solvent name.
- For ¹H, type *ss=8 nt=16* <rtn>
- Type *ga* <rtn> to start the experiment. If the concentration of the sample is low, you may want to increase nt.
- For ¹³C, use *nt=1e6 bs=8* <rtn>
- Type *ga* <rtn>. For 13C, after a few data blocks are complete (message: BS # completed), type *wft* to process. When sufficient S/N is obtained or your experiment time slot is finish, type *aa* to stop the experiment.
- To save, type *svf('Desire_name')* <rtn>

5. Reference Your Spectrum:

- Type *dscale* for to display ppm scale and find your solvent peak.
- With left mouse button (LC) click to the left of solvent peak. With right mouse (RC) button click to the right of solvent peak.
- Click **Expand**. Place cursor on top of solvent peak with **LC**.
- Type *nl rl(7.26p)* <rtn> for CDCl3.c

	¹ H ppm	¹³ C ppm	HOD ppm	Mel/Boil °C
CDCl ₃	7.24 (1)	77.23 (3)	1.5	-63.5/61
Acetone-d ₆	2.05 (5)	29.8 (7)	2.8*	-94/56.5
Benzene-d ₆	7.16 (1)	128.39 (3)	0.4	5.5/80.1
DMSO-d ₆	2.50 (5)	39.51 (7)	3.3*	3.81/101.42
CD ₃ CN	1.94 (5)	1.3 (7)	2.1*	-45/81.6
D ₂ O	4.80 (1)		4.8 (1)	3.81/101.42
CD ₃ OD	3.31 (5)	49.15 (7)		-97.8/64.7
CD ₂ Cl ₂	5.32 (3)	54.00 (5)	1.5	-95/39.75
DMF-d ₇	2.92 (5)	34.89 (7)		-61/153
	2.75 (5)	29.76 (7)		

6. Common VNMR Commands:

<i>aa</i> – abort acquisition	<i>aph</i> – auto-phase correction
<i>cexp(2)</i> – create experiment	<i>dscale</i> – display scale
<i>f</i> – display full spectrum	<i>full</i> – display spectrum in full
<i>ga</i> – acquire and process	<i>go</i> – acquire spectrum
<i>jexp1</i> – join experiment 1	<i>nl</i> – nearest line
<i>nt</i> – number of transients	<i>page</i> – send to plotter
<i>pap</i> – plot all parameters	<i>pl</i> – plot spectrum
<i>pscale</i> – plot scale	<i>res</i> – resolution of peak
<i>rl</i> – reference line: rl(7.24p)	<i>ss</i> – steady state scans
<i>su</i> – setup hardware param	<i>svf</i> – save FID
<i>svs</i> – save shims only	<i>unlock(2)</i> – unlock experiment 2
<i>vsadj</i> – vertical scale adjust	
<i>wft</i> – weighted Fourier transform	
<i>pl pap pscale page</i>	print plot with parameters
<i>plloo</i>	print plot with peak pick
<i>pl1d</i>	print plot with integration
<i>stackplot</i>	stack two plots together